

«Sticky Stuff» Adhesive Matting

The universal solution



- ▶ Universally applicable
- ▶ High shaking speeds
- ▶ Practical size and shape
- ▶ Temperature range 16–60°C
- ▶ Easy removal of flasks
- ▶ Easy to clean
- ▶ Can be disinfected
- ▶ Applicable to universal tables
- ▶ Easy to replace
- ▶ Unique – only from INFORS HT!



www.infors-ht.com

May be subject to technical amendments. Product images examples only.

INFORS HT

Unbeatable Quality

► Universally applicable

The INFORS HT adhesive matting «Sticky Stuff» is suitable for fastening of various types of culture vessels (maximum shaking speed depends on shape and bottom surface of vessel).

► High shaking speeds

Important for applications with high oxygen demand: up to 300 rpm are possible (depending on vessel).

► Practical size and shape

The convenient format of 200 x 200 mm makes placing small vessels and replacing the matting much easier.

► Temperature range 16–60°C

The «Sticky Stuff» shows very even properties over a wide temperature range.

► Easy removal of flasks

“Gently push/pull and wait a few seconds” – this way, culture vessels can be removed without effort. A very important safety aspect.

► Easy to clean

Simply brush off under running water – if required, use gentle dishwashing soap.

► Can be disinfected

Disinfection with quaternary ammonium compounds.

► Applicable to universal tables



The «Sticky Stuff» can be applied to existing universal trays with threaded holes.

► Easy to replace

Application without air bubbles is blindingly easy when the matting is wet.

► Unique – only from INFORS HT!

These positive features are only provided by INFORS HT green «Sticky Stuff».



Specification

Maximum shaking speeds: Only valid for flasks with 20% filling volume. Flasks and adhesive matting must be totally undamaged, clean, dry and free of grease. All information supplied without liability.

Schott Duran® glass Erlenmeyer flasks:

25–750 mL:	max. 200 rpm
1000 mL:	max. 250 rpm
2000 mL:	max. 250 rpm
3000 mL:	max. 300 rpm
5000 mL:	max. 250 rpm

Corning plastic Erlenmeyer flasks:

125–3000 mL:	max. 200 rpm
--------------	--------------

Corning plastic Fernbach flasks:

3000 mL:	max. 250 rpm
----------	--------------

Data valid for 50 mm shaking diameter. With 25 mm shaking diameter, slightly higher shaking speeds are possible.

Infors AG
Headoffice, Switzerland

Rittergasse 27
CH-4103 Bottmingen
T +41 (0)61 425 77 00
F +41 (0)61 425 77 01
info@infors-ht.com

For more information and your local sales office please visit:

www.infors-ht.com

INFORS HT