

Sound proof cabinets for ultrasonic baths

Sonation series USBB



**-30dB(A)
Over 90% reduction in noise**

Easy access

Space-saving

Robust construction

Active ventilation



- Noise reduction up to 30dB(A) which is a drop in noise by over 90%. The cavitation noise is practically eliminated.
- Available in three different sizes for mostly all ultrasonic baths.
- The cover can be completely opened, so the removal of the basket is very easy to handle.
- There is no additional space needed to open the cover, so the box can be placed and operated under the extractor hood or under lab cabinets.
- Control panel stays visible. Thanks to the transparent front, the progress of the washing program can be observed anytime.
- Equipped with fast access openings on the front and both sides to operate the ultrasonic bath or to drain the cleaning agent.
- Integrated fan to remove vapors.
- Recessed grips on the sides for easy transportation.
- Oil and water repellent and flame-retardant acoustic foam.
- Powder coated, robust and easy cleanable surface.
- Cushioned opening mechanism for smooth and easy handling.

Contact

Sonation GmbH
Alte Schulstr. 39
88400 Biberach
Tel.: +49 (0) 7351 30038-39
Fax.: +49 (0) 7351 30038-41

www.sonation.com

Models

The noise reduction boxes for ultrasonic baths are available in three sizes, S, M and L. The list below gives an overview over the supported ultrasonic baths. If your model is not listed, please contact Sonation and ask for the best fitting noise reduction box.

Bandelin

DT31
DT31H
DT52
DT52H
DT100
DT100H
DT102H
DT103H
DT106
DT156
DT156H
DT 255
DT 255 H
DT 510
DT 510 H
DT 512 H
DT 514
DT 514 H
DT 514 BH
DT 1028
DT 1028 H
DT 102 H-RC
DT 255 H-RC
DT 510 H-RC
DT 514 BH-RC
RK 31
RK 31 H
RK 52
RK52 H
RK100
RK100H
RK100 SH
RK102 H
RK103 H
RK106
RK156
RK156 BH
RK255
RK255 H
RK510
RK510 H
RK512 H
RK513
RK514

RK514 H
RK514 BH
RK1028
RK1028 H
RK102 CH
103 CH
RK255 CH
RK512 CH
Grand Instruments
XUB5
XUB12
XUB18
XUBA1
XUBA3

ELMA

S10/H
S15/H
S30/H
S40/H
S60/H
S100/H
S120/H
S180/H
P30 H
P60 H
P120 H
P180 H
Elmasonic Denta Pro
Elmasonic x-tra 30 H
Elmasonic x-tra 50 H
Elmasonic x-tra 70 H
Elmasonic x-tra 150 H
TI-H-5
T490 DH
T710 DH
T760 DH

Sonica

1200 M
2200 M
2200 M
2200 M
2200 M
2400 M

2400 MH
2400ETH
2400 EP
3200 M
3200 MH
3200 ETH
3200 EP
3200 L M
3200 L MH
3200 L ETH
3200 L EP
3300 L M
3300 MH
3300 ETH
3300 EP
4200 M
4200 MH
4200 ETH
4200 EP
4300 M
4300 MH
4300 ETH
4300 EP

Advantage Lab

AI04-01-230
AI04-03-230
AI04-06-230
AI04-12-230

Contact

Sonation GmbH
Alte Schulstr. 39
88400 Biberach
Tel.: +49 (0) 7351 30038-39
Fax.: +49 (0) 7351 30038-41